



## Morphine (MS Contin) Hapten/Protein Ratio=70:1 [BSA] (DAGF-171)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                     |                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------|
| Antigen Description | Morphine Bovine Serum Albumin (BSA) Conjugate                                           |
| Conjugate           | BSA                                                                                     |
| Format              | Liquid                                                                                  |
| Concentration       | Batch dependent - please inquire should you have specific requirements.                 |
| Size                | 1 mg                                                                                    |
| Buffer              | 10 mM Phosphate Buffered Saline, pH 7.4                                                 |
| Preservative        | 0.1% Sodium Azide                                                                       |
| Storage             | Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles. |

### BACKGROUND

|          |                                                                                                                                                                                |
|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Keywords | Chloride, Morphine;Contin, MS;Duramorph;Morphia;Morphine;Morphine Chloride;Morphine Sulfate;Morphine Sulfate (2:1), Anhydrous; Morphine Sulfate (2:1), Pentahydrate;MS Contin; |
|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|