



## Barbital/Barbitone [BSA] (DAGF-174)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Barbital Bovine Serum Albumin Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

### BACKGROUND

**Keywords** Barbital;Barbitone;Diemal;Diethylmalonylurea;Dormileno;Ethylbarbital;Medinal;Veronal