



## Buprenorphine [BSA] (DAGF-179)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                     |                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------|
| Antigen Description | Buprenorphine Bovine Serum Albumin Conjugate                                            |
| Conjugate           | BSA                                                                                     |
| Format              | Liquid                                                                                  |
| Concentration       | Batch dependent - please inquire should you have specific requirements.                 |
| Size                | 1 mg                                                                                    |
| Buffer              | 10 mM PBS, pH 7.4                                                                       |
| Preservative        | 0.1% Sodium Azide                                                                       |
| Storage             | Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles. |

### BACKGROUND

|          |                                                                                                                                                                           |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Keywords | 6029 M;6029-M;6029M;Buprenex;Buprenorphine;Buprenorphine Hydrochloride;Buprex;Hydrochloride, Buprenorphine;Prefin;RX 6029 M; RX-6029-M;RX6029M;Subutex;Temgésic;Temgesic; |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|