



Benzoylecgonine (COOH) Hapten/Protein Ratio=30:1 [BSA] (DAGF-183)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Benzoylecgonine-Bovine Serum Albumin Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Benzodiazepine;1,2-Benzodiazepine;benzodiazepin;1H-1,2-benzodiazepine;12794-10-4;benzodiazapine;UNII-M0Q7802G2B;D0CS6N;SCHEMBL8137;AC1L350D;CTK0I3497;DTXSID90155730;SVUOLADPCWQTTE-UHFFFAOYSA-N;M0Q7802G2B;Benzodiazepine (Aversion, anxiety);DB12537;HE031475;HE305111;LS-33991;20520-EP2269990A1
----------	---