



Cocaine (COC) [BSA] (DAGF-184)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Cocaine Bovine Serum Albumin Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Cocaine; Cocaine HCl; Cocaine Hydrochloride; HCl, Cocaine; Hydrochloride, Cocaine
----------	---