



Methamphetamine (MET) Hapten/Protein Ratio=50:1 [BSA] (DAGF-193)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Methamphetamine Bovine Serum Albumin Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords Deoxyephedrine;Desoxyephedrine;Desoxyn;Hydrochloride, Methamphetamine;Madrine;Metamfetamine;Methamphetamine;Methamphetamine Hydrochloride;Methylamphetamine;N Methylamphetamine;N-Methylamphetamine;