



Methamphetamine (MET) [BSA] (DAGF-194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|---|
| Antigen Description | Methamphetamine Bovine Serum Albumin Conjugate |
| Conjugate | BSA |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | PBS, pH 7.4 |
| Preservative | 0.02% Sodium Azide |
| Storage | Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles. |

BACKGROUND

| | |
|----------|--|
| Keywords | Deoxyephedrine;Desoxyephedrine;Desoxyn;Hydrochloride, Methamphetamine;Madrine;Metamfetamine;Methamphetamine;Methamphetamine Hydrochloride;Methylamphetamine;N Methylamphetamine;N-Methylamphetamine; |
|----------|--|