



Tricyclic Antidepressants (TCA) [BSA] (DAGF-199)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Tricyclic Antidepressant Bovine Serum Albumin (BSA) Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM PBS, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Chloralvan;Cohidrate;Dormal;Escr;HS;Hydral;Hydrate de chloral;Hydratedechloral;Hynos;Kessodrate;Kloralhydrat;Knockoutdrops;Lorinal;Lycoral;Tricyclic Antidepressant;TCA;TCAs
-----------------	--