



Pyrethroid [KLH] (DAGF-204K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Pyrethroids are axonic excitotoxins, the toxic effects of which are mediated through preventing the closure of the voltage-gated sodium channels in the axonal membranes. The sodium channel is a membrane protein with a hydrophilic interior.
--------------	---