



Recombinant HIV type 1 (MN) Glycoprotein 120 [FITC] (DAGH010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Baculovirus derived gp120 (HIV type 1 MN) conjugated to FITC.
Specificity	This fluoresce conjugated gp120 binds to murine monoclonal antibodies of defined epitope specificity and immobilized CD4 in Dot Blot assays.
Conjugate	FITC
Applications	Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence
Stability	At least 6 months at -75°C
Bio-activity	This fluoresce conjugated gp120 binds to CD4 bearing cells in FACS.
Format	Frozen Liquid
Size	50 µg
Preservative	None
Storage	-75 °C

BACKGROUND

Keywords	HIV type 1 gp120; Human Immunodeficiency Virus-1 gp120; HIV-1 gp120; Envelope glycoprotein GP120; gp120
----------	---