



## NFE2L2 (D77V mutant form) [His] (DAGH020)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	NFE2L2 protein is encoded by NFE2L2 gene. It is a basic leucine zipper protein that regulates the expression of antioxidant proteins that protect against oxidative damage triggered by injury and inflammation. Several drugs that stimulate the NFE2L2 pathway are being studied for treatment of diseases that are caused by oxidative stress.
<b>Species</b>	Human
<b>Purity</b>	>99% by SDS-PAGE
<b>Conjugate</b>	His
<b>Molecular Weight</b>	26 kDa
<b>Size</b>	100 µl
<b>Preservative</b>	None
<b>Storage</b>	-80°C

### BACKGROUND

<b>Keywords</b>	NFE2L2; nuclear factor, erythroid 2-like 2; NRF2; nuclear factor erythroid 2-related factor 2; HEBP1; NFE2-related factor 2; NF-E2-related factor 2; nuclear factor erythroid-derived 2-like 2; nuclear factor (erythroid-derived 2)-like 2; nuclear factor, erythroid derived 2, like 2
-----------------	--