



NFE2L2 (R34Q mutant form) [His] (DAGH021)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NFE2L2 protein is encoded by NFE2L2 gene. It is a basic leucine zipper protein that regulates the expression of antioxidant proteins that protect against oxidative damage triggered by injury and inflammation. Several drugs that stimulate the NFE2L2 pathway are being studied for treatment of diseases that are caused by oxidative stress.
Species	Human
Purity	>99% by SDS-PAGE
Conjugate	His
Molecular Weight	26 kDa
Size	100 µl
Preservative	None
Storage	-80°C

BACKGROUND

Keywords	NFE2L2; nuclear factor, erythroid 2-like 2; NRF2; nuclear factor erythroid 2-related factor 2; HEBP1; NFE2-related factor 2; NF-E2-related factor 2; nuclear factor erythroid-derived 2-like 2; nuclear factor (erythroid-derived 2)-like 2; nuclear factor, erythroid derived 2, like 2
-----------------	--