



Rat IgG2a kappa Isotype Control [PerCP/Cy5.5®] (DAGIC255)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	The SUL3869 immunoglobulin reacts with KLH. The isotype of this antibody is rat IgG2a, κ. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live and fixed mouse, rat and human tissues.
Immunogen	Trinitrophenol + KLH
Isotype	IgG2a, κ
Source/Host	Rat
Species Reactivity	N/A
Clone	SUL3869
Purification	Affinity chromatography
Conjugate	PerCP/Cy5.5
Applications	FCM
Preparation	The immunoglobulin was purified by affinity chromatography and conjugated with PerCP/Cy5.5 under optimal conditions. The solution is free of unconjugated PerCP/Cy5.5 and unconjugated immunoglobulin.
Format	Liquid
Concentration	0.2 mg/ml
Size	25 µg
Buffer	Phosphate-buffered solution, pH 7.2

Preservative	0.09% Sodium Azide
---------------------	--------------------

Storage	The antibody solution should be stored undiluted between 2°C and 8°C and protected from prolonged exposure to light. Do not freeze.
----------------	-------------------------------------------------------------------------------------------------------------------------------------

BACKGROUND

Introduction	<p>Isotype controls are primary antibodies that lack specificity to the target, but match the class and type of the primary antibody used in the application. Isotype controls are used as negative controls to help differentiate non-specific background signal from specific antibody signal. Depending upon the isotype of the primary antibody used for detection and the target cell types involved, background signal may be a significant issue in various experiments.</p>
---------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------
