



Vitamin D [Biotin] (DAGP-0201B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Synthetic Vitamin D antigen, Biotin conjugated
Purity	96.7% by HPLC (215 nm)
Conjugate	Biotin
Format	Purified, Lyophilized
Concentration	Not Applicable (Reconstitute in Ethanol)
Size	1 mg
Buffer	Lyophilized from an organic solvent.
Storage	Short time storage at -20 °C, long time storage at -80 °C. Shelf life 60 months