



Estriol-3 [HRP] (DAGP-2524H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Estriol(E3), also spelled oestriol, is a steroid, a weak estrogen, and a minor female sex hormone. It is one of three major endogenous estrogens, the others being estradiol and estrone. Levels of estriol in women who are not pregnant are almost undetectable.

Species N/A

Conjugate HRP

Applications ELISA

Format Liquid, Linking Position-3-CME position.

Size 1 ml

Buffer Supplied in Tris with Enzyme Stabilizers.

Preservative Anti-Bacterial Agent.

Storage Store at 2-8 °C.