



## Reverse triiodothyronine [BSA] (DAGP21819)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Reverse triiodothyronine (rT3) is an isomer of triiodothyronine (T3). It is the third-most common iodothyronine the thyroid gland releases into the bloodstream, the levels of rT3 will increase in conditions such as euthyroid sick syndrome because its clearance decreases while its production stays the same.
<b>Species</b>	N/A
<b>Conjugate</b>	BSA
<b>Applications</b>	LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.