



# 1-Aminohydantoin (AHD)-HRP Conjugate (DAGPY3201)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Antigen Description** 1-Aminohydantoin (AHD)-HRP conjugate is used for developing an ELISA assay for quantifying 1-Aminohydantoin (AHD). 1-Aminohydantoin (AHD)-HRP conjugate has minimal number of 1-Aminohydantoin (AHD) molecules on each HRP to maximize its assay sensitivity.

**Conjugate** HRP

**Format** Lyophilized powder

**Size** 1 mg

**Buffer** 10 mM PBS, pH 7.4

**Preservative** None

**Storage** Store at -20°C.

## BACKGROUND

**Keywords** 1-Aminoimidazolidine-2,4-dione; 1-Aminoimidazolidine-2,4-dione; 1-Amino hydantoin; 6301-02-6; Aminohydantoin