



Semicarbazide (SEM)-HRP Conjugate (DAGPY3205)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Semicarbazide (SEM)-HRP conjugate is used for developing an ELISA assay for quantifying Semicarbazide (SEM). Semicarbazide (SEM)-HRP conjugate has minimal number of Semicarbazide (SEM) molecules on each HRP to maximize its assay sensitivity.
Conjugate	HRP
Format	Lyophilized powder
Size	1 mg
Buffer	10 mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Keywords	aminourea; hydrazinecarboxamide; Carbamoylhydrazine; SEMICARBAZIDE, BASE; SEMICARBAZIDE; SEM
-----------------	--