



DHEA [BSA] (DAGPY6646)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DHEA is an important endogenous steroid hormone. It is the most abundant circulating steroid in humans, in whom it is produced in the adrenal glands, the gonads, and the brain, where it functions predominantly as a metabolic intermediate in the biosynthesis of the androgen and estrogen sex steroids.
Purity	> 95% pure
Conjugate	BSA
Applications	ELISA, Latex Immunoturbidimetric Assay
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	500 µg
Buffer	Supplied in liquid form dissolved in PBS pH 7.4
Preservative	None
Storage	Aliquot and store at -20 °C or -80 °C. Avoid repeated freeze/thaw cycles

BACKGROUND

Keywords	DHEA; Andrestenol; Androstenolone; Dehydroepiandrosterone; Dehydroisoandrosterone; Diandron; Diandrone; 5-dehydro-epiandrosterone; Gynodian; Prasterona; Prasterone; Prasteronum; Prestara; Psicosterone; Siscelar plus; trans-Dehydroandrosterone;
-----------------	---