



17B(3) ESTRADIOL [HRP] (DAGPYE2301)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	17beta-estradiol is the 17beta-isomer of estradiol. It has a role as an estrogen, a human metabolite, an EC 1.2.3.1 (aldehyde oxidase) inhibitor, a Daphnia magna metabolite and a mouse metabolite. It is a 17beta-hydroxy steroid and an estradiol.
Species	N/A
Conjugate	HRP
Applications	ELISA, Immunoassays
Format	Concentrate
Concentration	Please refer to the vial label for the specific concentration.
Size	500 µL
Buffer	Supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution.
Storage	Store at 2-8 °C for up to 3 months or -20 °C for long term storage
Ship	Wet ice