



Recombinant HCV Core, NS3, NS4, NS5 Antigen (Capture) (DAGRM008T)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS buffer, pH 7.4
Storage	Store < 2 weeks at 2-8°C. Long term store at -20°C. Aliquot to avoid multiple freeze/thaw cycles

BACKGROUND

Keywords HCV; Hepatitis C virus; Core antigen