



1-Aminohydantoin hydrochloride standard (DAGS028)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Use this reference material as calibrant for chromatography or any other analytical technique.
Purity	99.8%(Please refer to the certificate of analysis of each batch.)
Conjugate	N/A
Cas.No	2827-56-7
Applications	Reference material for laboratory use only.
Molecular Weight	151.55g/mol
Size	100 mg
Storage	The reference material should be stored in the original sealed container at the 20±4°C. After use the bottle should be tightly closed and protected from moisture.
Ship	Ship at room temperature.