



4-epi-Oxytetracycline standard (DAGS039)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This reference material is intended for use in a laboratory as a calibration and quality control standard or in method development for analytical techniques.
Conjugate	N/A
Cas.No	14206-58-7
Applications	Reference material for laboratory use only.
Molecular Weight	460.43g/mol
Format	Neat
Size	10 mg
Storage	The reference material should be stored in the original sealed container at the $-18\pm 4^{\circ}\text{C}$. After use the bottle should be tightly closed and protected from moisture.
Ship	Ship at room temperature.