



Azaperone standard (DAGS045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This Reference Material is intended for use in a laboratory as a calibration and quality control standard or in method development for analytical techniques.
Purity	98.2%(Please refer to the certificate of analysis of each batch.)
Conjugate	N/A
Cas.No	1649-18-9
Applications	Reference material for laboratory use only.
Molecular Weight	327.4g/mol
Format	Neat
Size	50 mg
Storage	The reference material should be stored in the original sealed container at the 20±4°C. After use the bottle should be tightly closed and protected from moisture.
Ship	Ship at room temperature.