



## Benzalkonium chloride standard (DAGS047)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                     |                                                                                                                                                                    |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Antigen Description | Use this reference material as calibrant for chromatography or any other analytical technique.                                                                     |
| Purity              | 94.6%(Please refer to the certificate of analysis of each batch.)                                                                                                  |
| Conjugate           | N/A                                                                                                                                                                |
| Cas.No              | 8001-54-5                                                                                                                                                          |
| Molecular Weight    | 227.77g/mol                                                                                                                                                        |
| Format              | Neat                                                                                                                                                               |
| Size                | 100 mg                                                                                                                                                             |
| Storage             | The reference material should be stored in the original sealed container at the 20±4°C. After use the bottle should be tightly closed and protected from moisture. |
| Ship                | Ship at room temperature.                                                                                                                                          |