



## Sulfapyridine standard (DAGS103)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                     |                                                                                                                                                                   |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Antigen Description | Use this reference material as calibrant for chromatography or any other analytical technique.                                                                    |
| Purity              | 99.58%(Please refer to the certificate of analysis of each batch.)                                                                                                |
| Conjugate           | N/A                                                                                                                                                               |
| Cas.No              | 144-83-2                                                                                                                                                          |
| Molecular Weight    | 249.29g/mol                                                                                                                                                       |
| Format              | Neat                                                                                                                                                              |
| Size                | 250 mg                                                                                                                                                            |
| Storage             | The reference material should be stored in the original sealed container at the 4±4°C. After use the bottle should be tightly closed and protected from moisture. |
| Ship                | Ship at room temperature.                                                                                                                                         |