



Magic™ Anti-KIM-1 monoclonal antibody, clone C2694N (DCAB-TJ218)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No cross-reactivity with TIM-3/HAVcr-2/KIM-3 or TIM-4/TIMD-4
Immunogen	Recombinant protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	C2694N
Purification	> 90% pure (SDS-PAGE). Protein A Chromatography
Conjugate	Unconjugated
Applications	Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended Suggested pair for testing (Capture - Detection): DCAB-TJ218 - DCAB-TJ219
Format	Purified, Liquid
Concentration	Lot Specific
Size	1 mg
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Preservative	0.095 % NaN ₃

Storage

Store at 2-8°C.

BACKGROUND

Introduction

Kidney injury molecule-1 (KIM-1) is a marker for renal proximal tubular damage, the hallmark of virtually all proteinuric, toxic and ischaemic kidney diseases.

Keywords

KIM-1; Kidney injury molecule-1

GENE INFORMATION

Entrez Gene ID

[26762](#)

UniProt ID

[Q96D42](#)
