



Anti-BCL2L1 monoclonal antibody, clone 3I23 (DCABH-10312)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to Bcl-XL
Antigen Description	Potent inhibitor of cell death. Inhibits activation of caspases (By similarity). Appears to regulate cell death by blocking the voltage-dependent anion channel (VDAC) by binding to it and preventing the release of the caspase activator, CYC1, from the mitochondrial membrane. Isoform Bcl-X(S) promotes apoptosis.
Specificity	This antibody recognizes a 27 kDa protein identified as Bcl-X protein. This antibody is highly specific to bcl-X and does not show cross reactivity with Bcl-2 alpha or bax protein.
Immunogen	Synthetic peptide corresponding to Human Bcl-XL aa 3-14. Sequence: C-QSNRELVVDFLS
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Mouse, Rat, Human, Pig
Clone	3I23
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr, IP, WB
Positive Control	Hodgkin's lymphoma.
Format	Liquid
Size	250 µl
Buffer	Preservative: 0.05% Sodium azide

Preservative	0.05% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	BCL2L1 BCL2-like 1 [Homo sapiens]
Official Symbol	BCL2L1
Synonyms	BCL2L1; BCL2-like 1; bcl-2-like protein 1; Bcl X; bcl xL; bcl xS; BCL2L; BCLX; PPP1R52; protein phosphatase 1; regulatory subunit 52; apoptosis regulator Bcl-X; protein phosphatase 1, regulatory subunit 52; BCLXL; BCLXS; Bcl-X; bcl-xL; bcl-xS; BCL-XL/S; D
Entrez Gene ID	598
Protein Refseq	NP_001182
UniProt ID	Q07817
Chromosome Location	20q11.21
Pathway	Amyotrophic lateral sclerosis (ALS), organism-specific biosystem; Amyotrophic lateral sclerosis (ALS), conserved biosystem; Apoptosis, organism-specific biosystem; Apoptosis, organism-specific biosystem; Apoptosis, conserved biosystem; Apoptosis, organism-specific biosystem; BH3-only proteins associate with and inactivate anti-apoptotic BCL-2 members, organism-specific biosystem;
Function	BH3 domain binding; identical protein binding; protein binding; protein heterodimerization activity; protein kinase binding;