



Anti-KLK9 monoclonal antibody (DCABH-12122)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KLK9 belongs to the kallikrein subgroup of serine proteases, which have diverse physiologic functions in many tissues. For background information on kallikreins, see MIM 147960.
Immunogen	A synthetic peptide of human KLK9 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	KLK9 kallikrein-related peptidase 9 [Homo sapiens]
Official Symbol	KLK9
Synonyms	KLK9; kallikrein-related peptidase 9; kallikrein 9; kallikrein-9; KLK L3; kallikrein-like protein 3;

KLKL3; KLK-L3; MGC50513;

Entrez Gene ID	284366
Protein Refseq	NP_036447
UniProt ID	Q2XQG6
Chromosome Location	19q13.33
Function	peptidase activity; serine-type endopeptidase activity;
