



Magic™ Anti-CDV Monoclonal antibody, Clone N29723DE2 (DCABY-4369)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CDV
Immunogen	Canine Distemper Virus antibody was raised in mouse using CDV as the immunogen.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	CDV
Clone	N29723DE2
Purification	Canine Distemper Virus antibody was purified by Protein A chromatography.
Conjugate	Unconjugated
Applications	ELISA Pr* Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4369 - DCABY-4370
Format	Liquid
Size	1 mg
Buffer	Supplied in PBS, pH 7.4, 0.1 % sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4 °C

BACKGROUND

Introduction

Canine distemper is a viral disease that affects animals in the families Canidae, Mustelidae, Mephitidae, Hyaenidae, Ailuridae, Procyonidae, Pinnipedia, some Viverridae and Felidae (though not domestic cats; feline distemper or panleukopenia is a different virus exclusive to cats). It is most commonly associated with domestic animals such as dogs and ferrets, although it can infect wild animals as well.

Keywords

Canine Distemper Virus; CDV;