



Magic™ Anti-E. coli O157 Monoclonal antibody, Clone N232035FD-G[FITC] (DCABY-4411)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	E. coli O157
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	E. coli
Clone	N232035FD-G
Purification	E. coli O157 antibody was purified by protein A chromatography
Conjugate	FITC
Applications	ELISA Pr*, IF, IHC, WB Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4445 - DCABY-4411
Format	Liquid
Size	100 µg
Buffer	Supplied in 0.01 M phosphate buffered saline buffer, pH 7.2 and contains 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at -20 °C until use

BACKGROUND

Introduction Serotype E. coli O157:H7 is a gram-negative bacterium that can produce a bloody diarrhea due

to toxins, especially Shiga (Vero) toxin, it secretes when it infects the human intestine.

Keywords

E. coli O157; E.coli; Escherichia Coli O157; E.coli O157 Escherichia Coli O157; E.coli O157;
