



Magic™ Anti-FABPs monoclonal antibody, clone N73956 (DCABY-4420)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FABP antibody was raised in mouse using human heart FABP as the immunogen
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	N73956
Purification	FABP antibody was purified by chromatography on protein A sepharose
Conjugate	Unconjugated
Applications	ELISA Pr*, Immunoassays Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4418 - DCABY-4420
Format	Liquid
Size	1 mg
Buffer	Supplied in PBS buffer, pH 7.4, with 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4 °C

BACKGROUND

Introduction FABP proteins are thought to facilitate the transfer of fatty acids between extra- and intracellular

membranes.

Keywords

Fatty acid-binding protein 3 Heart-type; Fatty acid-binding protein; H-FABP; Mammary-derived growth inhibitor; MDGI; Muscle fatty acid-binding protein; M-FABP; FABP3;

GENE INFORMATION

Entrez Gene ID

[2168](#)

UniProt ID

[P07148](#)
