



Magic™ Anti-FABPs monoclonal antibody, clone N80290 (DCABY-4421)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human
Immunogen	FABP antibody was raised in mouse using human heart FABP as the immunogen.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	N80290
Purification	FABP antibody was purified by chromatography on protein A Sepharose.
Conjugate	Unconjugated
Applications	ELISA Pr*, WB Suggested pair for testing: DCABY-4421 - DCABY-4422
Format	Liquid
Size	1 mg
Buffer	Protein A purified and supplied in PBS, pH 7.4, 0.1 % NaN ₃ .
Preservative	0.1% Sodium Azide
Storage	Store at 4 °C

BACKGROUND

Introduction

The fatty-acid-binding proteins (FABPs) are a family of carrier proteins for fatty acids and other lipophilic substances such as eicosanoids and retinoids. These proteins are thought to facilitate the transfer of fatty acids between extra- and intracellular membranes.

Keywords

FABP; Fatty Acid Binding Protein;

GENE INFORMATION

Entrez Gene ID

[2168](#)

UniProt ID

[P07148](#)