



# Magic™ Anti-Ficolin H monoclonal antibody, clone N92725-C [Biotin] (DCABY-4429)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Ficolin H antibody was raised in Mouse using Recombinant human H-ficolin as the immunogen
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	N92725-C
<b>Purification</b>	Ficolin H antibody was purified from cell culture supernatant by protein A/G chromatography
<b>Conjugate</b>	Biotin
<b>Applications</b>	ELISA Pr* Suggested pair for sandwich ELISA (Capture - Detection): <a href="#">DCABY-4428</a> - DCABY-4429
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	Supplied in 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
<b>Preservative</b>	15mM Sodium Azide
<b>Storage</b>	Store at 4-8 °C without exposure to light

## BACKGROUND

**Introduction** Ficolin H is a thermolabile beta-2-macroglycoprotein found in all human serum and is a member

of the ficolin/opsonin p35 lectin family.

---

**Keywords**

Ficolin H; FCN3; H-Ficolin; Hakata antigen; Ficolin-3; Ficolin 3;

---

## GENE INFORMATION

**Entrez Gene ID**

[8547](#)

---

**UniProt ID**

[O75636](#)

---