



Mouse Anti-Human GLP-1(aa 1-36) monoclonal antibody, clone N92690 (DCABY-4439)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Anti-GLP-1(1-36)amide (N-terminal-extension specific) monoclonal antibody
Specificity	DCABY-4439 is specific for the N-terminal region (residues 1-6) of the GLP-1. No cross-reactitvity with GLP-1 (7-36) amide or (7-37). DCABY-4439 reacts with the low activity enzyme GLP-1(1-36)amide, but not with either of the activated forms GLP-1(7-36)amide and GLP-1(7-37)amide
Immunogen	Synthetic peptide
Isotype	lgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	N92690
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Epitope	Residues 1-6 of GLP-1
Format	Liquid
Concentration	Lot specific
Buffer	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Preservative	15mM Sodium Azide
Storage	4-8°C without exposure to light. No precautions necessary during handling.
Ship	Wet ice

BACKGROUND

Introduction

Glucagon-like peptide 1(7-36)amide (GLP-1(7-36)amide) is the principal active form of GLP-1, the other being GLP-1(7-37). GLP-1 is a peptide hormone of the glucagon family, produced by the L cells of the intestinal mucosa from the same prohormone as glucagon. The active forms are potent stimulators of glucose-dependent insulin secretion. The sequence of GLP-1 is fully conserved in all mammalian species examined so far.

 $\begin{tabular}{ll} \textbf{Keywords} & GCG; glucagon; GLP1; GLP2; GRPP; preproglucagon \\ \end{tabular}$

GENE INFORMATION

 Entrez Gene ID
 2641

 UniProt ID
 P01275