



Mouse Anti-Human GLP-1(aa 7-36) monoclonal antibody, clone N92637 [Biotin] (DCABY-4443)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Anti-GLP-1 (GLP-1(7-36)amide, C-terminal specific) monoclonal antibody
Specificity	DCABY-4443 is specific for the amidated C-terminus of the peptide and does not react with GLP-1(7-37)
Immunogen	Synthetic GLP-1(7-36)amide coupled to carrier
Isotype	lgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	N92637
Purification	Protein A purified
Conjugate	Biotin
Applications	ELISA(Det)
Epitope	C-terminal epitope of GLP-1(7-36)amide
Format	Liquid
Concentration	Lot specific
Buffer	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl
Preservative	15mM Sodium Azide

45-1 Ramsey Road, Shirley, NY 11967, USA

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Storage	4-8°C without exposure to light. No precautions necessary during handling.
---------	--

Ship Wet ice

BACKGROUND

Introduction Glucagon-like peptide 1(7-36)amide (GLP-1(7-36)amide) is the principal active form of GLP-1,

> the other being GLP-1(7-37). GLP-1 is a peptide hormone of the glucagon family, produced by the L cells of the intestinal mucosa from the same prohormone as glucagon. The active forms are potent stimulators of glucose-dependent insulin secretion. The sequence of GLP-1 is fully

conserved in all mammalian species examined so far.

Keywords GCG;glucagon;GLP1;GLP2;GRPP;preproglucagon

GENE INFORMATION

Entrez Gene ID	<u>2641</u>
UniProt ID	<u>P01275</u>

Email: info@creative-diagnostics.com