



Magic™ Anti-LCN2 monoclonal antibody, clone N92645 (DCABY-4574)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NGAL is a component of granules in neutrophils from tissues that are normally exposed to microorganisms and is upregulated during inflammation.
Immunogen	NGAL antibody was raised in Mouse using Recombinant monkey NGAL as the immunogen
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Monkey
Clone	N92645
Purification	NGAL antibody was purified from cell culture supernatant by protein A/G chromatography
Conjugate	Unconjugated
Applications	ELISA Pr*, IHC Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4574 - DCABY-4583
Format	Liquid
Size	200 µg
Buffer	Supplied in 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Preservative	15mM Sodium Azide
Storage	Store at 4-8 °C without exposure to light

GENE INFORMATION

Gene Name	LCN2 lipocalin 2 [Macaca mulatta]
Official Symbol	LCN2
Synonyms	LCN2; lipocalin 2;
Entrez Gene ID	697208
UniProt ID	H9H4D4
Chromosome Location	chromosome: 15