



Magic™ Anti-Staphylococcus type D Staphylococcus Enterotoxin D Monoclonal antibody, Clone N230635 (DCABY-4637)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to Staphylococcal enterotoxin type D
Target	Staphylococcus type D Staphylococcus Enterotoxin D
Immunogen	Staphylococcus Enterotoxin Type D antibody was raised in mouse using Staphylococcal enterotoxin D as the immunogen.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Staphylococcus
Clone	N230635
Purification	Staphylococcus Enterotoxin Type D antibody was purified by Protein A or G affinity purification
Conjugate	Unconjugated
Applications	ELISA Pr* Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4636 - DCABY-4637
Format	Liquid
Size	500 µg
Buffer	Supplied as a liquid in PBS, pH 7.4, 3mM KCl, 1.5mM KH ₂ PO ₄ , 140mM NaCl, 8.0mM Na ₂ HPO ₄

Preservative	None
Storage	Store at 4 °C for short term or -20 °C for longer storage. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction	<p>An enterotoxin (not to be confused with endotoxin) is a protein toxin released by a microorganism in the intestine. Enterotoxins are chromosomally encoded exotoxins that are produced and secreted from several bacterial organisms. They are often heat-stable, and are of low molecular weight and water-soluble. Enterotoxins are frequently cytotoxic and kill cells by altering the apical membrane permeability of the mucosal (epithelial) cells of the intestinal wall. They are mostly pore-forming toxins (mostly chloride pores), secreted by bacteria, that assemble to form pores in cell membranes. This causes the cells to die.</p>
Keywords	Staphylococcus Enterotoxin Type D; Enterotoxin Type D; Enterotoxin D;