



Magic™ Anti-Serpinf2 monoclonal antibody, clone N60867 (DCABY-4665)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Alpha 2 antiplasmin antibody was raised in mouse using purified alpha-2 antiplasmin as the immunogen.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	N60867
Conjugate	Unconjugated
Applications	ELISA Pr* Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4665 - DCABY-4667
Format	Liquid
Size	500 µg
Buffer	Supplied in lyophilized form in PBS, pH 7.4.
Preservative	None
Storage	Store at 4 °C until reconstitution. Following reconstitution aliquot and freeze at -20 °C for long term storage. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction	Alpha 2-antiplasmin (or alpha 2-antiplasmin or plasmin inhibitor) is a serine protease inhibitor
---------------------	--

(serpin) responsible for inactivating plasmin, an important enzyme that participates in fibrinolysis and degradation of various other proteins. This protein is encoded by the SERPINF2 gene.

Keywords	alpha 2 Antiplasmin; alpha-2 Antiplasmin; A2AP;
-----------------	---

GENE INFORMATION

Entrez Gene ID	5345
-----------------------	----------------------

UniProt ID	P08697
-------------------	------------------------
