



Magic™ Anti-Serpinf2 monoclonal antibody, clone N60869 (DCABY-4667)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Alpha 2 antiplasmin antibody was raised in mouse using purified alpha-2 antiplasmin as the immunogen.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	N60869
Conjugate	Unconjugated
Applications	ELISA Pr* Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4665 - DCABY-4667
Format	Liquid
Size	500 µg
Buffer	Supplied in lyophilized form from PBS, pH 7.4, with 0.02% NaN3 and 20mg/ml mannitol.
Preservative	0.02% Sodium Azide
Storage	Store at 4 °C until reconstitution. Following reconstitution aliquot and freeze at -20 °C for long term storage. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction Alpha 2-antiplasmin (or alpha 2-antiplasmin or plasmin inhibitor) is a serine protease inhibitor

(serpin) responsible for inactivating plasmin, an important enzyme that participates in fibrinolysis and degradation of various other proteins. This protein is encoded by the SERPINF2 gene.

Keywords	alpha 2 Antiplasmin; -; alpha-2 Antiplasmin; A2AP;
-----------------	--

GENE INFORMATION

Entrez Gene ID	5345
-----------------------	----------------------

UniProt ID	P08697
-------------------	------------------------
