



# Magic™ Anti-N. gonorrhoeae Monoclonal antibody, Clone N2800OH3 (DCABY-4773)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Gonorrhea
Target	N. gonorrhoeae
Immunogen	Neisseria gonorrhoeae antibody was raised in Mouse using NGH as the immunogen
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N. gonorrhoeae
Clone	N2800OH3
Conjugate	Unconjugated
Applications	ELISA Pr*, Lateral Flow Pr* Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4773 - <a href="#">DCABY-4772</a>
Format	Liquid
Size	1 mg
Buffer	Supplied in 10 mM PBS, 50 mM NaCl, 0.05% NaN <sub>3</sub>
Preservative	0.05% Sodium Azide
Storage	Store at 2-8 °C for short term storage. For long term storage aliquot and store at -20 °C, avoid repeated freeze/thaw cycles

## BACKGROUND

**Introduction**

*Neisseria gonorrhoeae*, also known as gonococci (plural), or gonococcus (singular), is a species of Gram-negative coffee bean-shaped diplococci bacteria responsible for the sexually transmitted infection gonorrhea.

---

**Keywords**

*Neisseria gonorrhoeae*; Gonorrhea;

---