



## User's Manual

# Human Transcription factor 20 ELISA kit



DEIA-BJ1176



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Human Transcription factor 20 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Transcription factor 20. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

### Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

### Detection Range

1.0-25 ng/mL

### Sensitivity

0.1 ng/mL