



User's Manual

Rat A Disintegrin And Metalloprotease 8 ELISA Kit



DEIA-BJ2039



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Rat A Disintegrin And Metalloprotease 8 ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the A Disintegrin And Metalloprotease 8. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

50-1000 pg/mL

Sensitivity

0.1 ng/mL