



## User's Manual

# Rat Plasma levels of terminal complement complex ELISA Kit

REF

DEIA-BJ2053



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Rat Plasma levels of terminal complement complex ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Plasma levels of terminal complement complex. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

### Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

### Detection Range

25-500 ng/mL

### Sensitivity

1.0 ng/mL