



User's Manual

Rat Cholinergic Receptor, Muscarinic 2 ELISA Kit



DEIA-BJ2080



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Rat Cholinergic Receptor, Muscarinic 2 ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Cholinergic Receptor, Muscarinic 2. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

0.5 -10 ng/mL

Sensitivity

0.1 ng/mL