



User's Manual

Mouse Dihydropyrimidinase Like Protein 3 ELISA Kit

REF

DEIA-BJ2423



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA



Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com



Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Mouse Dihydropyrimidinase Like Protein 3 ELISA Kit kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Dihydropyrimidinase Like Protein 3. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

2.5-50 ng/mL

Sensitivity

0.1 ng/mL