



**User's Manual**

# Mouse Growth Associated Protein 43 ELISA Kit

REF

DEIA-BJ2466



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

Mouse Growth Associated Protein 43 ELISA Kit kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Growth Associated Protein 43. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

### Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

### Detection Range

0.5 -10 ng/mL

### Sensitivity

0.1 ng/mL