



User's Manual

Mouse Lctin like oxLDL Receptor 1 ELISA Kit



DEIA-BJ2505



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Mouse Lctin like oxLDL Receptor 1 ELISA Kit kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Lctin like oxLDL Receptor 1. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

100-2500 pg/mL

Sensitivity

1.0 pg/mL