



**User's Manual**

# **Mouse Platelet Derived Growth Factor Receptor Like Protein ELISA Kit**

**REF**

**DEIA-BJ2540**



**96T**

**RUO**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

Mouse Platelet Derived Growth Factor Receptor Like Protein ELISA Kit kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Platelet Derived Growth Factor Receptor Like Protein. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

### Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

### Detection Range

250-5000 pg/mL

### Sensitivity

1.0 pg/mL