



User's Manual

Rabbit Basic Fibroblast growth factor ELISA Kit



DEIA-BJ2657



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Rabbit Basic Fibroblast growth factor ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Basic Fibroblast growth factor. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

0-1000 pg/mL

Sensitivity

1.0 pg/mL